



# CATALOGUE



## FULL EUROPEAN CERTIFICATION

Manufacturer: Balerina-Czech, s.r.o.  
Přestavlcká 1447, Mnichovo Hradiště  
Czechia



   [balerina.cz](https://www.balerina.cz)



# MEDICAL FFP2 RESPIRATOR — PFHM731


**FILTRATION EFFICIENCY IS 99% FOR THE PARTICLES OF THE SIZES OF 120 NANOMETRES**  
(the size of virus SARS-CoV-2)

Complies also with the legal requirements on medical face masks (the requirements on breathability, microbial cleanliness, liquid resistance, filtration efficiency) and biocompatibility (the product is not cytotoxic, is not skin irritant, does not show any sensibilization potential).


The product is in conformity with **Regulation (EU) 2016/425** on personal protective equipment and also in conformity with **Regulation (EU) 2017/745** on medical devices.

Parameter	Legal Requirement	Balerina PFHM731
Filtration Efficiency	min. 94 %	99,6 %
Breathing resistance (walk)	max. 70 Pa	40,2 Pa
Breathing resistance (run)	max. 240 Pa	116,5 Pa
Exhalation resistance	max. 300 Pa	176,7 Pa


# SIZES — PFHM731



S




S/M



M/L


# FFP2 RESPIRATOR — PFHM731

DESIGN 5




M/L

DESIGN 6



S S/M M/L

DESIGN 8



S S/M M/L

DESIGN 14




S S/M

DESIGN 15



S S/M M/L


DESIGN 22



M/L


# FFP2 RESPIRATOR — PFHM731

DESIGN 28



M/L

DESIGN 32




M/L

DESIGN 35



M/L

DESIGN 36




S S/M M/L

DESIGN 37



M/L

DESIGN 39



M/L

# FFP2 RESPIRATOR — PFHM731

DESIGN 43




S S/M M/L

DESIGN 44




S S/M M/L

DESIGN 47



S S/M

DESIGN 48




M/L

DESIGN 49



S S/M M/L


DESIGN 51



S S/M M/L

# FFP2 RESPIRATOR — PFHM731

DESIGN 53




S S/M M/L

DESIGN 55




M/L

DESIGN 56




S S/M M/L

DESIGN 57



M/L

DESIGN 58



M/L

DESIGN 62



S S/M M/L

# FFP2 RESPIRATOR — PFHM731

DESIGN 63



S S/M M/L

DESIGN 76




M/L

DESIGN 77




S S/M M/L

DESIGN 87



M/L

DESIGN 89



M/L

DESIGN 90




M/L




# FFP2 RESPIRATOR — PFHM731

DESIGN 92




M/L

DESIGN 96




M/L

DESIGN 100




M/L

DESIGN 101




M/L

DESIGN 102



M/L


DESIGN 105



M/L

# FFP2 RESPIRATOR — PFHM731

DESIGN 106



M/L

DESIGN 108



S S/M M/L

DESIGN 110




M/L

DESIGN 116




M/L

DESIGN 118



M/L


DESIGN 119



M/L

# FFP2 RESPIRATOR — PFHM731

**DESIGN 121**




**S S/M M/L**

**DESIGN 122**




**M/L**

**DESIGN 125**



**M/L**

**DESIGN 150**




**S S/M M/L**

**DESIGN 171**



**S S/M M/L**


**DESIGN 183**



**S S/M**

# FFP2 RESPIRATOR — PFHM731

DESIGN 185




S S/M

DESIGN 186




S S/M

DESIGN 187




S S/M M/L

DESIGN 198




S S/M M/L

DESIGN 200



M/L

DESIGN 203



M/L


# FFP2 RESPIRATOR — PFHM731

DESIGN 221




M/L

DESIGN 224




M/L

DESIGN 227




M/L

DESIGN 229




M/L

DESIGN 231



M/L

DESIGN 234



S S/M M/L

# FFP2 RESPIRATOR — PFHM731

**DESIGN 238**



**S S/M M/L**

**DESIGN 247**



**S S/M M/L**

**DESIGN 256**




**S S/M**

**DESIGN 263**




**S S/M**

**DESIGN 266**



**S S/M M/L**


**DESIGN 268**



**M/L**

# FFP2 RESPIRATOR — PFHM731

**DESIGN 274**




**M/L**

**DESIGN 314**




**M/L**

**DESIGN 332**



**M/L**

**DESIGN 999**



**M/L**



# MEDICAL FFP2 RESPIRATOR — PFHM931

**FILTRATION EFFICIENCY IS 99% FOR THE PARTICLES OF THE SIZES OF 120 NANOMETRES**  
(the size of virus SARS-CoV-2)

Complies also with the legal requirements on medical face masks (the requirements on breathability, microbial cleanliness, liquid resistance, filtration efficiency) and biocompatibility (the product is not cytotoxic, is not skin irritant, does not show any sensibilization potential).


The product is in conformity with **Regulation (EU) 2016/425** on personal protective equipment and also in conformity with **Regulation (EU) 2017/745** on medical devices.

Parameter	Legal Requirement	Balerina PFHM931
Filtration Efficiency	min. 94 %	<b>99,8 %</b>
Breathing resistance (walk)	max. 70 Pa	<b>47 Pa</b>
Breathing resistance (run)	max. 240 Pa	<b>154 Pa</b>
Exhalation resistance	max. 300 Pa	<b>262 Pa</b>




# FFP2 RESPIRATOR — PFHM931

DESIGN 8




DESIGN 28




DESIGN 35




DESIGN 53



DESIGN 56




DESIGN 57




# FFP2 RESPIRATOR — PFHM931


DESIGN 87




DESIGN 89



DESIGN 100




DESIGN 101



DESIGN 105



DESIGN 106



# FFP2 RESPIRATOR — PFHM931


**DESIGN 110**




**DESIGN 116**



**DESIGN 118**




**DESIGN 125**



**DESIGN 198**




**DESIGN 203**




# FFP2 RESPIRATOR — PFHM931


DESIGN 224




DESIGN 227



DESIGN 231



DESIGN 234



DESIGN 238



DESIGN 314



# FFP2 RESPIRATOR — PFHM931

**DESIGN 332**



**DESIGN 999**





# MEDICAL FFP3 RESPIRATOR — PFHM933

**FILTRATION EFFICIENCY IS 99% FOR THE PARTICLES OF THE SIZES OF 120 NANOMETRES**  
(the size of virus SARS-CoV-2)


Complies also with the legal requirements on medical face masks (the requirements on breathability, microbial cleanliness, liquid resistance, filtration efficiency) and biocompatibility (the product is not cytotoxic, is not skin irritant, does not show any sensibilization potential).

The product is in conformity with **Regulation (EU) 2016/425** on personal protective equipment and also in conformity with **Regulation (EU) 2017/745** on medical devices.


Parameter	Legal Requirement	Balerina PFHM933
Filtration Efficiency	min. 99 %	99,6 %
Breathing resistance (walk)	max. 100 Pa	41 Pa
Breathing resistance (run)	max. 300 Pa	135 Pa
Exhalation resistance	max. 300 Pa	228 Pa

# FFP3 RESPIRATOR — PFHM933


DESIGN 5



DESIGN 6




DESIGN 8



DESIGN 15



DESIGN 28




DESIGN 48




# FFP3 RESPIRATOR — PFHM933


DESIGN 53




DESIGN 56




DESIGN 87




DESIGN 89



DESIGN 101




DESIGN 105





# FFP3 RESPIRATOR — PFHM933

DESIGN 106




DESIGN 110




DESIGN 116



DESIGN 118



DESIGN 125




DESIGN 198




# FFP3 RESPIRATOR — PFHM933


**DESIGN 203**



**DESIGN 224**




**DESIGN 231**




**DESIGN 314**



**DESIGN 332**



**DESIGN 999**





# MEDICAL FFP2 RESPIRATOR — PFHM722


Complies also with the legal requirements on medical face masks (the requirements on breathability, microbial cleanliness, liquid resistance, filtration efficiency) and biocompatibility (the product is not cytotoxic, is not skin irritant, does not show any sensibilization potential).

The product is in conformity with **Regulation (EU) 2016/425** on personal protective equipment and also in conformity with **Regulation (EU) 2017/745** on medical devices.

Parameter	Legal Requirement	Balerina PFHM722
Filtration Efficiency	min. 94 %	99,6 %
Breathing resistance (walk)	max. 70 Pa	23,5 Pa
Breathing resistance (run)	max. 240 Pa	87,2 Pa
Exhalation resistance	max. 300 Pa	150 Pa

# FFP2 RESPIRATOR — PFHM722

**DESIGN 314**




**M/L**

**DESIGN 332**



**M/L**

**DESIGN 999**



**M/L**